

# ARCHITECT SYSTEM

## ASSAY PARAMETERS (Preliminary\*)

## $\beta$ -Hydroxybutyrate LiquiColor®

Serum/Plasma SI Units

GENERAL PARAMETERS		
Name:	BHOB	
Assay#:	#	
Assay Version:	1	
Assay Type:	Photometric	
Assay Availability:	Enabled	
REACTION DEFINITION		
Reaction Mode:	END UP	
Primary Wavelength:	500	
Secondary Wavelength:	-	
Last Read Required:	33	
Absorbance Range:	0.000-0.000	
Sample Blank Type:	Self	
Blank Assay:	-	
Main Read Time:	15 – 31	
Flex Read Time:	-	
Blank Read Time:	14 - 16	
REAGENT / SAMPLE		
Reagent:	BHOB	
R1 Reagent Volume:	215	
R1 Water Volume:	0	
R1 Dispense Mode:	Type 0	
R2 Reagent Volume:	36	
R2 Water Volume:	0	
R2 Dispense Mode:	Type 0	
Diluent Dispense Mode:	Type 0	
Standard Sample Volume:	6.0	1:100
1:3 Sample Volume:	2.0	1:2.95
1:10 Sample Volume:	20.0	1:10.00
Diluted Sample Volume :	6.0	
Diluent Volume:	180	
	*	
VALIDITY CHECKS		
Reaction C heck Type:	None	
Reaction Time A Range:	-	
Calculation Limit:	-	
Read Time B Range	-	
Minimum Absorbance:	-	
Rate Linearity %:	-	
CALIBRATION PARAMETERS		
Calibration Method:	Linear	
Use Cal From:	-	
Full Interval Hours:	72	
Adjustment Type:	None	
Expected Cal Factor:	0.00	
Exp. Cal Factor Tolerance %	0	
Span:	Blank	
Max. Curve Fit:	0.00	
Calibrator Set Name	BHOB	
Factor:		
Adjustment Interval Hours:	0	
Default Ordering Type:	Full	
	0.0000-	
Blank Absorbance Range:	0.0000	
Span Absorbance Range:	0.00 – 0.00	
Cal Level (Water) Conc.:	0	
Cal Level (Water) Vol.:	6.0	
Cal Level (BHOB1) Conc.:	1.00	
Cal Level (BHOB1) Vol.:	6.0	
RESULT PARAMETER		
Normal Range:	0.02 – 0.27	mmol/L
Result Decimal Places:	2	
Correlation Factor::	1.0000	
Intercept:	0.0000	

### Order information

**Cat. No. 2440-058**

### Notes

- Please refer to the package insert of the  $\beta$ -Hydroxybutyrate LiquiColor® for detailed information about the test on the following:

Clinical Relevance  
 Method and Principle  
 Composition and Stability of the Reagents  
 Specimens  
 Calibrators and Controls  
 Performance Characteristics regarding  
   - Measuring Range  
   - Specificity/Interferences  
   - Sensitivity/Limit of Detection  
   - Precision (Reproducibility, Repeatability)  
   - Method Comparison

Reference Ranges  
 Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by  
 Stanbio Laboratory  
 1261 North Main Street • Boerne, Texas 78006 USA  
 www.stanbio.com

**\*Customers should be aware that this protocol has been generated by a user of our assay on an Abbott Architect analyzer and has not been validated by Stanbio Laboratory or the instrument manufacturer.**

Architect is a registered trademark of Abbott Diagnostics.

- #) Data entry by the user  
 \*) Enter calibration or standard value and position  
 \*\*) Factor to be checked by a calibration serum