

**STANBIO REAGENT APPLICATION (VALIDATED)****Procedure No. 2440****B-HYDROXYBUTYRATE LIQUICOLOR<sup>®</sup>****COBAS MIRA CHEMISTRY SYSTEMS\*****\*Roche Diagnostic Systems****Nutley, NJ****PLEASE REFER TO THE PACKAGE INSERT FOR ADDITIONAL INFORMATION REGARDING THIS REAGENT.****Package Sizes Available**

Catalog No. 2440-058 (1 x 58mL)

**Reagent Preparation/Stability**Reagents are supplied ready to use. **DO NOT USE CLEANER IN APPLICATION.**Reagents are stable stored at 2-8°C until expiration date on their respective labeling. **Once opened, contamination must be avoided.****Linearity**

When perform as directed, this method is linear to 4.5 mmol/L.

**MIRA INSTRUMENT SETTINGS****GENERAL**

Measurement Mode: Absorb. (1)  
Reaction Mode: R – S-SR1 (1)  
Calibration Mode: Slope Avg. (1)  
Reagent Blank: Reag/Dil (2)  
Cleaner: No (1)  
Wavelength: 500 nm (2)  
Decimal Position: 2  
Unit: mmol/l (21)

**ANALYSIS**

Sample Dil.: Name: H2O (00)  
Post. Dil.: Factor: No (Space)  
Conc. Factor: No (Space)  
Sample: Cycle: (1)  
Volume: uL (3.0)  
Diluent: uL (20.0)  
Reagent: Cycle: (1)  
Volume: uL (105)  
Start R1 Cycle: 2  
Volume: uL (18)  
Dilution Name: H2O  
Volume: uL (60)

**CALCULATION**

Sample Limit: No (Space)  
Reaction Direction: Increase (1)  
Check: On (1)  
Factor: 1.000  
Offset: 0.000  
Test Range: Low (0)  
High (6)  
Normal Range (37°C): Low (.9)  
High: (6)  
Number of Steps: 1  
Calc. Step A: Endpoint (3)  
Readings: First: 1  
Last: 12  
Reaction: Limit: No

**CALIBRATION**

Calib. Interval: On Request (3)  
Reagent Range: Low: No  
High: No  
Blank Range: Low: No  
High: No  
Calibrator: Cup Position (#)  
Cal-1: +  
Replicate: Triplicate (3)  
Deviation: % (10)

**CONTROL**

CS1 Pos/Low/Assign/High (#)  
CS2 Pos/Low/Assign/High (#)  
CS3 Pos/Low/Assign/High (#)

# User Defined

\* The Factor must be checked by a calibrator serum or control

**Expected Values**

0.02 to 0.27 mmol/l

Due to geographic or instrument variation, it is recommended that each laboratory establish the normal range for its population.

**MIRA S INSTRUMENT SETTINGS****GENERAL**

Measurement Mode: Absorb (1)  
Reaction Mode: R – S-SR1 (1)  
Calibration Mode: Slope Avg. (1)  
Reagent Blank: Reag/Dil (2)  
Cleaner: No (1)  
Wavelength: 500 nm (2)  
Decimal Position: 2  
Unit: mmol/l (21)

**ANALYSIS**

Post Dil.: Factor: No (Space)  
Conc. Factor: No (Space)  
Sample: Cycle: (1)  
Volume: uL (3.0)  
Dilution: Name: H2O (00)  
Volume: uL (20.0)  
Reagent: Cycle: (1)  
Volume: uL (105)  
Start R1 Cycle: 2  
Volume: uL (18)  
Dilution Name: H2O  
Volume: uL (60)

**CALCULATION**

Sample Limit: No (Space)  
Reaction Direction: Increase (1)  
Check: On (1)  
Factor: 1.000  
Offset: 0.000  
Test Range: Low (0)  
High (6)  
Normal Range (37°C): Low (.9)  
High: (6)  
Number of Steps: 1  
Calc. Step A: Endpoint (3)  
Readings: First: 1  
Last: 12  
Reaction: Limit: No

**CALIBRATION**

Calib. Interval: On Request (3)  
Reagent Range: Low: No  
High: No  
Blank Range: Low: No  
High: No  
Standard: Position (#)  
Std-1: +  
Std-2: No (Space)  
Std.-3: No (Space)  
Replicate: Triplicate (3)  
Deviation: % (10)

**CONTROL**

CS1 Pos/Low/Assign/High (#)  
CS2 Pos/Low/Assign/High (#)  
CS3 Pos/Low/Assign/High (#)

# User Defined

\* The Factor must be checked by a calibrator serum or control

**FOR TECHNICAL ASSISTANCE CALL TOLL-FREE 1-800-531-5535 OR (830) 249-0772**