

# Stanbio Laboratory

## $\beta$ -Hydroxybutyrate LiquiColor®

# ORTHO 5,1/5600 VITROS

### Reagent Order information

Ref. No. Kit size  
2440-058 R1 1 x 50 mL + R2 1 x 8.5 mL  
1 x 3 mL Std.

### Controls/Linearity Order information

TDM/ $\beta$ -Hydroxybutyrate Tri-level Controls, Ref. No. 2460-605  
TDM/ $\beta$ -Hydroxybutyrate Bi-level Controls, Ref. No. 2465-605  
TDM/ $\beta$ -Hydroxybutyrate Linearity Standards, Ref. No. 2450-604

### Method

$\beta$ -Hydroxybutyrate

### Reagent preparation and stability

The reagents are ready-to-use and stable up to the end of the indicated month of expiry, if contamination is avoided and stored at 2 – 8 °C. On-board stability have not been determined.

### Specimen

Serum or plasma collected with EDTA, heparin or sodium fluoride. Avoid hemolysis!  
Levels are stable at least 7 days at 2 - 8 °C.

### Components in the test

**R1:**  $\beta$ -hydroxybutyrate dehydrogenase  
Diaphorase

**R2:** NAD  
INT  
Oxalate

**Std:** 1.0 mmol/L

### Notes

- See interfering substances in package insert.
- The reagents contain sodium azide (0.095 %) as preservative. Do not swallow! Avoid contact with skin and mucous membranes!

### Normal Range

0.02 – 0.27 mmol/L

### Instrument setting

Full Assay Name: B-Hydroxybutyrate Date:  
Short Assay Name: BHOB Fluid Type: Serum  
Assay Model Type: End Point Template: EPT1 R1-s-R2  
Cal Model Type: Linear Calibrator bottle: 2  
Reagent Reps per Cal: 3

Standard Dilution Factor: 1.0 Diluent: Water  
Reflex Dilution Information  
Reflex Dilution: On Dilution Factor: 2.0 Reduction Factor: 1.0

Samples Indices Check: Enabled  
Threshold Limits  
Hemolysis: 1000 Icterus: 25 Turbidity: 800

Reporting Type: Quantitative Units: mmol/L  
Significant Digits: 3 Precision Digits: 2  
Slope: 1.00 Intercept: 0.00  
CuveTip Expiration Time: 35 Temperature Sensitive: No

Ranges  
Reference: 0.02 – 0.27  
Supplementary: 0.00 – 900000000  
Reportable: 0.00 – 4.50

Initial Abs. Limits: -0.200 – 2.700 Blank Abs. Limit: -0.200 - 2.700

Monotonicity: Increase  
Max Response High: 3.000 Max Response Low: -3.000  
Min Response High: 3.000 Min Response Low: -3.000

	Reportable Conc.	Triple Read Limit
Reportable Min:	0.00	400
Critical Conc.:	2.25	8.00
Reportable Max:	4.50	8.00

Kit Lot	Bottle Number	Dilution Factor	Cal. Value	Response Range
	1	1.0	0.00	0.20000
	2	1.0	1.00	0.20000

Protocol Steps

1. Reagent	Volume (uL) = 143.3, Pack/Bottle=UD01/A
2. Incubation	Seconds = 190.00
3. Sample	Volume (uL) = 4.0
4. Incubation	Seconds = 23.75
5. Read	Wavelength = 510nm
6. Incubation	Seconds = 4.75
7. Reagent	Volume (uL) = 24.0, Pack/Bottle=UD01/B
8. Incubation	Seconds = 313.50
9. Read	Wavelength = 510nm

### Note: UDA Cartridges:

Volume A = 14.7 mL

Volume B = 2.8 mL

This will yield approximately 80 tests.