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Stanbio RaPET® ASO Procedure No. 1125

Latex agglutination slide test for the Qualitative and Semi-quantitative determination of Antistreptolysin-O in non-diluted serum.

Summary and Principle

In infections caused by B-hemolytic streptococci, Streptolysin-O is liberated from the bacteria stimulating production of Antistreptolysin-O antibodies (ASO). The presence and level of these antibodies in a serum specimen may reflect the nature and severity of infection.¹

The Stanbio RaPET® ASO uses a stabilized buffered suspension of polystyrene latex particles that have been coated with Streptolysin-O.² When the latex reagent is mixed with a serum specimen containing antibodies to Streptolysin-O, agglutination occurs. The latex reagent has been produced so that agglutination will take place only when the level of antibodies to Streptolysin-O is at/or greater than 200 IU/mL, a level determined to be indicative of disease by epidemiological and clinical studies.³

Reagents

ASO Latex Reagent, Ref. No. 1126 (White Cap) Suspension of polystyrene latex particles coated with Streptolysin-O in a saline buffer.

ASO Positive Control, Ref. No. 1127 (Red Cap) Stabilized human serum, containing more than 200 IU/mL of Antistreptolysin-O.

Negative Control, Ref. No. 1192 (Green Cap) Stabilized human serum, containing less than 100 IU/mL of Antistreptolysin-O.

Glycine-Saline Buffer(20X) Concentrate, Ref. No. 1191 Solution of glycine and sodium chloride, pH 8.2 ± 0.1.

Warnings and Precautions: For In Vitro Diagnostic Use.

The controls used in this kit have been tested by an FDA-approved method and found non-reactive for the presence of HBsAg and antibody to HIV. While these methods are highly accurate, no test can offer complete assurance that infectious agents are absent. This material, as well as all patient samples, should be handled as though capable of transmitting infectious disease. The United States Food and Drug Administration recommends such samples be handled at Center for Disease Control's Biosafety Level 2.

Reagents in this kit contain sodium azide as a preservative. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.

Reagent Preparation and Stability

Reagents are stable when stored at 2-8°C until the expiration date shown on their respective labels. Prepare Glycine-Saline Buffer by adding the Concentrate (5mL) to a 100mL volumetric flask and diluting to the 100mL mark with distilled or deionized water.

Do Not Freeze!

Materials Provided

Glass Slide - 6 cell Disposable
Pipette/Mixers

Material Required (Not Provided)

Timer Rotator (optional)

Light Source

Test Tubes & Rack (Semi-quantitative test only) Serological Pipettes (Semi-quantitative test only)

Specimen Collection and Preparation

It is recommended that serum only be used. Samples must have clotted completely and contain no particulates nor traces of fibrin after thawing. Do

not heat inactivate test sera or controls. Avoid repeated freeze/thawing of specimens.

Sample Stability: Serum specimens are reportedly stable up to 8 days at 2-8°C and for up to 3 months if they are frozen (at/or below -25°C) within 24 hours of venipuncture and are not repeatedly thawed and re-frozen.

Interfering Substances: As in all serological test, hemolytic, lipemic or turbid sera may cause incorrect results and should not be used. Use only a clean, dry slide washed in mild detergent and rinsed thoroughly with distilled water.

Limitations

Reaction times greater than 2 minutes can lead to false positive results (due to drying effect). Very lipemic sera can also cause non-specific reactions.

Strength of agglutination is not indicative of the ASO concentration in the sample. In the qualitative test procedure, weak reactions may occur with markedly elevated concentrations. In cases of greatly increased ASO titer (more than 2000 IU/mL), agglutination may be inhibited because of antibody excess (prozone effect). When such high ASO concentrations are to be expected, the sample should be tested diluted. The results of this test should always be interpreted in the light of clinical and other laboratory findings. Negative results do not rule out the diagnosis of acute rheumatic fever or poststreptococcal glomerulonephritis.

Procedure

A. Qualitative Method


1. Bring all reagents and patient samples to room temperature.
2. Gently shake the reagent vial to disperse and suspend the latex particles.
Do Not Shake Vigorously!
3. Using a dispensing pipette provided, add one drop of the non-diluted patient sample(s) onto a separate cell(s) on the glass slide. Also place one (1) drop of the positive and negative control on separate cells of the test slide.
4. Gently mix the contents of the latex reagent including the contents of the glass dropper. Fill the dropper with the well-mixed latex suspension and place one (1) drop next to the drop of serum sample(s) on each of the separate cells in use.
5. Mix both drops with the mixer provided, covering the whole surface of each of the cell(s).
6. Tilt the slide back and forth for 2 minutes manually so that the mixture rotates slowly inside the cells or place the slide on a automatic rotator set at 80 – 100 rpm.
7. After the 2 minutes, examine each cell(s) for the absence or presence of agglutination.

Quality Control

A positive control and a negative control should be included in each test series. A positive control will produce a coarse agglutination against a clear background, while a negative control will produce a smooth, homogeneous suspension.

B. Semi-Quantitative Method

1. This reaction can be used to estimate the ASO concentration using a dilution series. The patient sample is diluted with the Glycine-Saline Buffer as shown:

Thereafter, each is tested with the latex as described under section "A. Qualitative Method."

Results

A. Qualitative Method

Agglutination identifies a ASO concentration greater than 200 IU/mL ± 20% in the sample.

Those samples that do not show agglutination contain ASO concentrations less than 200 IU/mL.

B. Semi-Quantitative Test Results

The highest sample dilution which still shows distinct agglutination is reported. The ASO content of the patient sample is taken from the following table:

Dilution of Agglutination	Concentration of ASO-IU/mL (± 20%)
1:2	≥ 400
1:4	≥ 800
1:8	≥ 1600
1:16	≥ 3200

Expected Values⁴

A detectable level of 200IU/mL Antistreptolysin-O antibodies is usually regarded as the normal upper limit, since less than 15-20% of healthy individuals demonstrate titers higher than 200IU/mL when their sera are assayed. Elevated ASO titers may be associated with ankylosing spondylitis, glomerulonephritis, scarlet fever and tonsillitis. Increased ASO levels are generally not found in sera of patients with Rheumatoid Arthritis except during acute episodes.

Performance Characteristics⁵

A comparison of 354 serum specimens were tested comparing this method with another commercially available procedure. Equivalent results were obtained on the qualitative method on 349 of the 354 random specimens. The overall agreement is approximately 99%.

References

1. Rantz L. D., DiCaprio JM, Randall E. Am. J. Med. Sci., 24, 1952.
2. Halbert, SP. Ann. N.Y. Acad. Sci., 103, 1027:1051; 1963.
3. Klein GL, Applied Microbiology, 21:999, 1971.
4. Klein GC: Manual of Clinical Immunology ASM, 264-273: 1976.
5. Stanbio Laboratory Data

STANBIO LABORATORY, LP DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES OF THE MERCHANTABILITY AND FITNESS PERTAINING TO THIS PRODUCT WHICH ARE NOT EXPRESSLY DETAILED IN THIS PACKAGING INFORMATION OR A WRITTEN AGREEMENT BETWEEN THE BUYER AND SELLER OF THIS PRODUCT.

STANBIO LABORATORY, LP MAINTAINS THAT THIS PRODUCT CONFORMS TO THE INFORMATION CONTAINED IN THIS INSERT. PURCHASER MUST DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS PARTICULAR USE. USE ONLY IN ACCORDANCE WITH LABELING INSTRUCTIONS.

Index of Symbols			
	Attention, see instructions for use		Tests per kit
	For <i>in vitro</i> diagnostic use only		Use by
	Store between temperature indicated		Lot No.
	European Representative		Flammable
	Manufactured by		Do not reuse
	Reference No.		Contents

Stanbio Laboratory
1261 North Main Street
Boerne, Texas 78006, USA
Toll Free: (800)-531-5535
Phone: (830) 249-0772
Fax: (830) 249-0851
e-mail: stanbio@stanbio.com
www.stanbio.com

EC REP
BDS International Diagnostics GmbH
Anselm-Feuerbach-Str. 12,
68723 Schwetzingen, Germany
e-mail: info@bds-international.de