

DiaSpect Tm



DiaSpect Tm

Immediate hemoglobin results
from a hand-held analyzer

- ◆ Results in less than two seconds
- ◆ Reagent-free cuvettes
- ◆ Internal battery for 40 days
or 10,000 tests continuous use
- ◆ 2 ½ years cuvette shelf life

Fast, easy to use and reliable in all climate conditions

Fast and accurate	• Results in less than two seconds • Measuring range: 1.2 - 25.5 g/dL • Precision (within run): CV ≤ 1% • Factory calibrated against HiCN reference method according to ICSH • Photometric method with broad spectrum and compensation for turbidity and scattering • No additional calibration required
Easy to use	• 'Always on' technology. The analyzer is always ready for use • Automatic self-check between every measurement • Long cuvette shelf life: 2 ½ years from date of production • Cuvettes unaffected by humidity or temperature and can be used until expiration date after opening the bag
Practical and portable	• Sample volume: 10 µL • Capillary or venous blood (EDTA and lithium heparin tubes can be used) • Maintenance free • Operating ambient temperature: 10 to 35°C (50 to 90°F) • Storage temperature: 0 to 50°C (32 to 122°F) • Transport temperature: -30 to 70°C (-22 to 158°F) up to a maximum of 24 hours. • Integrated rechargeable lithium-ion battery (40 days or 10,000 tests of continuous use) • Compact: 150 (d) x 90 (w) x 40mm (h) / 5.9 (d) x 3.5 (w) x 1.5 (h) inches • Lightweight: 180 g / 0.4 lb • Connect to PC via USB 2.0



1 Collect blood sample



2 Insert the filled cuvette and press down gently



3 Results appear in less than two seconds

U.S.A. Sales

EKF Diagnostics
1261 North Main Street
Boerne, TX 78006

☎ 1 800 531 5535
✉ sales2@ekfdiagnostics.com
www.ekfusa.com

Manufacturer

EKF-diagnostic GmbH
Ebendorfer Chaussee 3
39179 Barleben
Germany

EKF | Diagnostics
for life