



## Economical, accurate and consistent urine chemistry testing

- ◆ Minimal training time makes testing simple
- ◆ Data storage function recalls the last 2,000 test results
- ◆ Up to 120 tests per hour
- ◆ Lockout functions
- ◆ Built-in thermal printer

# State of the art electronics to accurately and efficiently read the Uri-Chek® 10SG urinalysis reagent strips.

Easy to use	<ul style="list-style-type: none"> <li>• Minimal training time</li> <li>• Recalls the last 2,000 test results</li> <li>• Up to 120 tests per hour in continuous test mode</li> <li>• QC and user lockout programmable</li> </ul>
Practical & portable	<ul style="list-style-type: none"> <li>• Built-in thermal printer</li> <li>• Barcode reader available</li> </ul>
Accurate & reliable	<ul style="list-style-type: none"> <li>• Reflectance Photometer Methodology</li> <li>• Photosensitive diode detection</li> <li>• 60 test/hr (single test mode)</li> <li>• 120 test/hr (continuous mode)</li> </ul>
Technical	<ul style="list-style-type: none"> <li>• Strip incubation: 1 min</li> <li>• Detection wavelengths: 525 nm and 635 nm</li> <li>• Power Source: 100-240V AC, 50/60Hz, 35W</li> <li>• Weight: 5.73 lbs / 2.6 kg</li> <li>• Dimensions (LxWxH): 10.7" x 10.6" x 5.7" (27cm x 27cm x 14cm)</li> </ul>

## Ordering information

Reference No.	Description
G1200-001	<b>Uri-Trak® 120 Urine Analyzer</b> Includes Analyzer , Strip Holder, Printer Paper Rolls (2), Fuses, (2.0A)(2), Power Cord, Quick Guide
1280-100	<b>Uri-Chek® 10SG Urinalysis Reagents Strips</b> Parameters: Urobilinogen, Glucose, Ketone, Bilirubin, Protein, Nitrite, pH, Blood, Specific Gravity and Leukocytes. 100 strips

Distributed by

Stanbio  
 1261 North Main Street  
 Boerne, Texas 78006  
 USA Toll-Free: 1 866 782 6246  
 Tel: +1 830 249 0772

[www.ekfusa.com](http://www.ekfusa.com)

**EKF** | Diagnostics  
 for life