

## Uri-Trak® 120M Urine Analyzer



### Uri-Trak® 120M

High accuracy and reliability  
for multiple parameters\*

- Automatic calibration before each test
- Up to 120 tests per hour in Continuous Mode
- Minimal training required
- 510(k) Cleared and CLIA Waived
- Data Management ready

\* Urobilinogen, Glucose, Ketones, Bilirubin, Protein, Nitrite, pH, Blood, Specific Gravity, and Leukocytes



**EKF**  
Diagnostics


User friendly operation	<ul style="list-style-type: none"> <li>• Semi-automatic operation with calibration before each test</li> <li>• 15 different test combinations per analyzer</li> <li>• Multiple language options: English, French, Italian, Portuguese and Spanish</li> <li>• High throughput: <ul style="list-style-type: none"> <li>◦ Single Test Option: 60 tests/hour</li> <li>◦ Continuous Test Option: 120 tests/hour</li> </ul> </li> <li>• 510(k) Cleared and CLIA Waived</li> </ul>
Accurate results with well-established and proven technology	<ul style="list-style-type: none"> <li>• Uses Reflectance photometry with photosensitive diode detection at wavelengths of 525 nm and 635 nm</li> <li>• Just 1 minute for strip incubation time</li> <li>• Parameters available: <ul style="list-style-type: none"> <li>◦ Glucose, Bilirubin, Ketone (Acetoacetic acid), Specific Gravity, pH, Blood, Protein, Urobilinogen, Leukocytes, Ascorbic Acid, Nitrite</li> </ul> </li> </ul>
Multiple data management options	<ul style="list-style-type: none"> <li>• Built in thermal printer or can connect and print to external printer</li> <li>• Data can transfer using either USB or RS232C options</li> <li>• Operator ID and Patient ID can be input either through manual entry or the optional RS232C Barcode Reader (up to 20 characters)</li> <li>• Memory can store up to 2000 results</li> </ul>
Small footprint	<ul style="list-style-type: none"> <li>• Compact 275 mm (l) x 183 mm (w) x 139 mm (h) (10.8" x 7.2" x 5.5")</li> <li>• Lightweight 1.85 kg (4.1 lbs) without batteries or power supply</li> </ul>

Catalog No.	Name
G1300-001	Uri-Trak® 120M Analyzer
1280-100	Uri-Chek® 10SG Urinalysis Reagent Strips
G1340-001	Data Transfer Cable
G1250-001	Barcode Reader
G1260-001	Printer, Thermal Paper
G1270-001	Printer, Sticker Paper

Distributed by

**Manufacturer**  
Stanbio Laboratory  
1261 North Main Street  
Boerne, Texas 78006

 (800) 531-5535  
 sales2@ekfdiagnostics.com

 [www.ekfdiagnostics.com](http://www.ekfdiagnostics.com)